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This quick guide does not substitute the careful reading of the user manual. For complete information about the use of this device, including safety warnings, charging the battery and the correct measurement for the monitored parameters, see the user manual found on the CD that accompanies the equipment.

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English

Warnings

IMPORTANT: this device must only be operated by qualified technical personnel. Before using, read the manual attentively.

WARNING: for defibrillation, cardioversion, pacing and AED modes, we do not recommend the use in patients under 01 (one) year old. For monitoring modes such as ECG, SpO₂, NIBP and EtCO₂, its use is possible with suitable accessories to these patients.

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ATTENTION: The patient's movements can confuse the correct detection of the rhythm and delay the therapy. Do not perform any movement with the patient and keep him a stopped position during rhythm analysis.

WARNING: DualMax can be used by patients over 01 (one) year old, regardless of their weight.

WARNING: risk of explosion if the equipment is operated in the presence of flammable liquids or gases.

ELECTRICAL SHOCK HAZARD: NEVER OPEN THE DEVICE. Each and every repair must be performed by Instramed's authorized technical centers.

WARNING: the patient must be placed on non-conductive surfaces, do not use wet or metallic surfaces and, if necessary, dry the chest before applying the shock.

WARNING: do not touch the patient, the equipment, the accessories nor any metallic or conductive surface which is in contact with the patient during the defibrillation.

WARNING: the patient must be completely still during the cardiac rhythm analysis phase (aed mode), do not give cardiac massage at this point.

ATTENTION: do not use DualMax or its accessories in the presence of MRI equipment.

This equipment was projected to offer resistance to electromagnetic interferences. However, the functioning of this device can be affected in the present of strong sources of electromagnetic-interference or radio-frequency, such as mobile phones, communicator radios, etc.

If the precision of measurements seems to be incorrect, first check the vital signs of the patient and then check the functioning of the DualMax.

WARNING: always check the general state of the equipment, the battery and the accessories before using it.

Before installing the equipment verify if there are any abnormalities or damage caused by mishandling during transportation.

NOTICE: the DualMax must only be used as a complement to assess the patient's physiological conditions. It must be accompanied by constant analysis of the patient's clinical status and symptoms.

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WARNING: The use of the DualMax is restricted to one patient at a time.

NOTICE: the applied parts (electrodes, sensors, armbands, etc.) are protected against defibrillation discharge; during discharge there may be baseline variation.

WARNING: when the DualMax is operated in monitor mode, it can be used with other electromedical equipment simultaneously connected to the patient, provided that the other equipment are in compliance with the safety standards.

WARNING: the conductive parts of the electrodes and connectors associated with the applied parts, including the neutral electrode, must not come into contact with other conductive parts, including the ground wire.

NOTICE: avoid connecting the patient to several items of equipment at the same time. The limits of current leakage may be exceeded.

NOTICE: the applied parts intended to come into contact with the patient have been evaluated and comply with the directives and principles of ISO 10993-1.

NOTICE: when removing the equipment from its package, carefully verify if there is any abnormality or visible damage in the device or its accessories, caused by impact or mishandling during transportation. In case of irregularities, please contact Instramed.

NOTICE: disposable accessories and any other components must be disposed of according to the norms of hospital waste disposal.

WARNING: Do not modify this equipment without authorization from Instramed.

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DualMax does not have support system and its carrying handle must not be used for this purpose.

This equipment does not have reminder signals.

Precautions, restrictions and warnings

The DualMax is a device built according to NBR and IEC standards and therefore offers total safety for patient and operator. However, all safety precautions described below must be followed.

 The monitor's operation can be affected by the presence of electromagnetic power sources, such as electrosurgical equipment and computer tomography (CT).

1 - ECG

1 - To guarantee protection against the effects of a defibrillation, use only the patient cable that accompanies the equipment.

2 - If the monitor is used simultaneously with an electric scalpel, position the ECG electrodes as far as possible from the RF current route, between the surgical field and the neutral card. Do not use needle type ECG electrode during surgical procedures.

2 - SpO₂

1 - The operation of this device may be affected by the presence of electromagnetic energy sources, such as electrosurgical or CT equipment. It also may be damaged by the presence of strong ambient light. If necessary, protect the area of the sensor with a surgical towel.

2 - Any dyes injected into the blood stream, such as methylene blue, indocaine green, indigo carmine and fluorescein, may affect the SpO₂ reading precision. The presence of dyshemoglobin, such as carboxyhemoglobin (in consequence of carbon monoxide poisoning) or methemoglobin (in consequence of sulfonamide's treatment) may affect the SpO₂ reading precision.

Adverse or side effects

Superficial burns may occur on the patient's skin in the area in contact with the electrodes. To minimize the effect of the disposable paddles, apply them directly after removal from the protection envelope and attach them firmly to the patient's skin.

The skin must be dry, or electric current leakage may occur, increasing the burn's area and reducing the efficiency of the treatment.

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Possibility of reddish skin at the defibrillation electrodes exposition place, due to the high voltage applied.

Possibility of superficial skin burns. To minimize the effect, in the case of adhesive pads, apply them immediately after removal of the protective envelope and securely attach to the patient's skin. The patient's skin must be dry.

Possibility of reduction of treatment efficiency. The patient's skin must be dry, otherwise the electric discharge may leak. Do not apply conductive gel outside the area of contact with the pads to avoid electrical leakage.

Do not apply conductive gel outside the area of contact with the pads to avoid electrical leakage.

Possibility of erythematous skin at the pacemaker electrodes exposition place, especially during prolonged use.

Possibility of discomfort and pain, with the use of high levels of pacemaker current. In case of high levels of current, it is recommended that the patient is anesthetized or unconscious.

Possibility of eczematous skin, due to the use of non-biocompatible electrodes. Always use accessories with technical recommendation described in this User Guide.

Possibility of reddish and/or bruised skin at the application place (thorax) by the use of CPR MAESTRO. It is recommended for cases of resuscitation maneuvers of long duration, the use of a gauze between the skin and the CPR MAESTRO.

Possible DISCOMFORT and REDUCED BLOOD CIRCULATION in NIBP measuring. The operator should always ensure that the applied pressure or the quantity of measures did not effect in the blood circulation or discomfort to the user. In this case, the quantity of measures should be decreased and the arm should be changed.

DISCOMFORT or COMPLICACIONES DUE TO THE PRODUCT SMELL in EtCO₂ monitoring. Some patients have been with nausea or headaches due to the cannula smell or discomfort with the cannula positioning. In this cases, it is reversible with just changing the cannula or readjusting the positioning.

Device care

- Do not place the equipment where it may fall on the patient. Do not lift the equipment by its cables or connections.
- Place cables connected to the patient in order to restrict the possibility of strangulation.
- Keep the defibrillator in a dry environment, avoiding places that allow liquids to spill over the monitor. Do not use the defibrillator if it is wet or excessively humid.
- Always keep the equipment and its accessories clean and well maintained.

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- If you suspect a fall or external damage, do not use the equipment.

IMPORTANT RECOMMENDATIONS

- Never immerse in liquid and never spill liquid of any kind on any part of the equipment.
- Do not use any other cleaning products not recommended by this manual.

NEVER sterilize any parts of the equipment, regardless of the sterilization method, as this would damage the mechanical structure and compromise the product's operation.

The recommendations above will ensure that the device will withstand, without damage or deterioration of safety factors, the necessary cleaning process.

Grounding

-  GROUNDING IS ESSENTIAL TO PROTECT THE OPERATOR AND PATIENT AGAINST ELECTRICAL DISCHARGE ACCIDENTS. IN THE ABSENCE OF ADEQUATE GROUNDING, DANGEROUS CURRENTS MAY CIRCULATE FROM THE EQUIPMENT BOX IF THERE IS AN INTERNAL ELECTRICAL DEFECT. GROUNDING MUST BE PERFORMED ACCORDING TO ABNT NORMS FOR ELECTRICAL INSTALLATIONS (NBR 13534/1995).

Electromagnetic compatibility

Warning:

Installing the DualMax requires special precautions concerning electromagnetic compatibility according to the information contained in this manual.

Mobile and portable RF communications equipment, such as mobile phones, can affect the DualMax's functioning.

Portable RF communication equipment (including peripherals such as cables and external antennas) should not be used within 30 cm of any part of DualMax, including cables specified by Instramed. Failure to do so may cause performance degradation.

Maximum length of accessories cables - In compliance with electromagnetic compatibility requirements:

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Warning:

Using cables, transducers and accessories different from the ones specified above, except for the ones sold by Instramed as replacement pieces, may result in emission increase or immunity decrease of the equipment.

The DualMax must not be used too close to or piled over other equipment.

When the equipment is used in a surgical procedure simultaneously with an electric scalpel, there is a risk of burnouts if a defect in the connection of the neutral electrode of the high frequency equipment matches a defect in the DualMax's ECG socket. This type of accident will only occur when the defects occur simultaneously, as the DualMax's ECG socket is electrically protected against risks of burnouts, being completely insulated.

The actions to be taken to prevent adverse events to the patient and the operator due to electromagnetic disturbances, during the useful life are:

- Ensure minimum distance, according to the table on page
- Electrical immunity - General, of an RF emitting source (radio frequency).
- The cables and accessories must also keep this distance.
- Do not use this product in conjunction with electrical scalpel.
- Do not use this product in conjunction with MRI devices.

The essential performance of the DualMax is understood as the ANALYSIS of the patient's heart rate, the DEFIBRILLATION, CARDIOVERSION and PACEMAKER applied to the patient, and the monitoring of the PATIENT

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in electrocardiogram, oximetry, non-invasive pressure, respiration, and invasive pressure. The performance of the DualMax is designed and verified for the absence of an unacceptable risk.

If performance is lost or degraded due to electromagnetic disturbances

- The ECG signal may be interfered with and cardiac rhythm analysis may be compromised as long as electromagnetic disturbances persist.
- The SPO₂ curve can be disturbed, besides the possibility of alteration in the measured value.
- Non-invasive pressure measurement may suffer alterations on the measured values.
- The Invasive pressure curve may suffer interferences.

Qualified users

Shall be considered qualified users those who have a degree in Medicine.

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When Non-Invasive Blood Pressure (NIBP) is present

- 1 Adult armband (cuff)*.
- 1 armband (cuff) extensor*.

When Invasive Blood Pressure is present

- 1 Utah Medical invasive blood pressure accessories kit*, consisting of:
 - 1 Organizer.
 - 1 Clamp.
 - 1 Electrical cable.
 - 1 Sampling line.
 - 1 Transducer*.

When Capnography (EtCO₂) is present

- 1 Oridion capnography accessories kit*, consisting of:
 - 1 Adult tracheal sampling line.
 - 1 Child tracheal sampling line.
 - 1 Adult nasal sampling line.

When PRINTER is present

- 1 White thermosensitive paper bobbin 58 mm x 15 m x 40 mm.

* Accessory with separate CE certificate.

OBS: the items already included may be acquired for replacement as extra items.

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 English

ATTENTION: observe the battery charge maintenance instructions.

First use

Before using the DualMax for the first time, the equipment must receive a full battery charge. In order to do this, the equipment needs to be connected to an electric current for at least eight hours.

Occasional use

Even when disconnected (stand-by), the DualMax executes internal routines checking the status of the equipment. Despite of this procedure entailing a low power consumption, the battery charge may be consumed.

Therefore, whenever the device has not been connected to an electric current for more than 20 days, it is advisable to execute a full battery charge. If this procedure is not performed, there is a risk of draining the battery and consequently being unable to use the DualMax in its portable configuration (not connected to the electric current).

Storage

The battery must be removed from the equipment in case it is stored or not used.

Replacement

Every battery has a determined lifetime, which is the possible quantity of full charge and discharge cycles, without loss of performance. When the appliance has a drop in performance of the battery, with low autonomy, request a new unit from Instramed technical assistance.

The battery can be replaced following the procedures described on the chapter "Care and maintenance".

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 **ATTENTION:** the equipment is capable of operating with only one battery connected.

 **ATTENTION:** when servicing the batteries, it is obligatory to change only one battery at a time, leaving the other battery in the device, so that it can continue to operate at that moment.

 More information on this topic is available in the DualMax User Manual present on the CD which accompanies the product.

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1. Touch screen.
2. Indication of AC Mains, battery charging and product operation status (see item "Indication of AC Mains, battery charge and product status" forward, in this chapter).
3. Print: quick access button (see item "Quick access buttons", forward, in this chapter).
4. E-Jog control: access to the equipment's general settings and menu navigation.
5. Quick access buttons (see item "Quick access buttons" forward, in this chapter).

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6. **Disarm:** disarms the stored charge.
7. **Sync:** activates synchronization.
8. **Selector key:** turns the equipment on and off; selects the operation mode.
9. **Charge:** charges the electrodes with the selected charge.
10. **Shock:** applies shock.
11. **Alarm LEDs:** indicate HIGH PRIORITY ALARM (RED) OR MEDIUM PRIORITY (YELLOW).

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Screen

The LCD screen displays graphic and numeric information used in ECG and SpO₂, defibrillation and others.

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
The e-Jog control is used to access all of the DualMax's functions, such as set alarms, define information displayed on the screen, alter parameters, etc.

ROTATE: rotating allows the user to select or change information and navigate all menus. It operates similarly to a computer mouse.

PRESS: works similarly to the enter button on a computer, confirming the selection.

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
Scale from 1 to 360 J: Allows the user to select the desired energy charge.

Monitor mode: used to monitor ECG, SpO₂, NIBP, EtCO₂, and RESP parameters, as in a multiparametric monitor.

Pacemaker mode (*): enables external pacemaker.

Position : Turns off the equipment.

AED position(*): enables external automatic defibrillator mode.

Auto Seq mode: enables charge auto-sequencing.

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NOTE: the equipment does not defibrillate in pacemaker and monitor modes. The pacemaker will only work in pacemaker mode.

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OP (*) Check your equipment's configuration. This item is optional and may not be present in all commercialized equipment.

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Quick access buttons

Fast lead change: enables quick access to change ECG leads.

Freeze: Freezes graphic screen signals for a closer examination (same Cancel button).

OP **Print (when available):** press once to print a quick report. For continuous printing, simply

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press the button for 3 seconds.

Pause audio: press the button quickly to deactivate ALL sound alarms for 60 seconds.

NIBP (when available): starts or suspends the functionality of the Non-Invasive Blood Pressure Measurement.

When the NIBP (optional) parameter is not present in the device, this button has no function.

Power and battery charging indications

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
- Power connected:** when the LED is on, it indicates that the equipment is connected to a power source or an external battery.
- Battery charging:** when the LED is on, it indicates that the battery is charging.
- Indicates that the equipment is ready to be used.

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
OBS: when the equipment is connected to an electric current, the LEDs will light up indicating the beginning of charging, even if the DualMax is inoperative.

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1. **Printer***: printer for thermosensitive paper. It prints electrocardiograms and events.
2. **Pads**: the pads accompanying the DualMax must be placed on top of the equipment and must be properly connected to the adult adapter.
3. **Transport handle**: Handle used to carry the equipment.
4. **Temperature connector***: Temperature connector YSI 400 standard (channel 1).
5. **Temperature connector***: Temperature connector YSI 400 standard (channel 2).
6. **ECG connector**: Connector for ECG cables. Depending on the parameters present in the equipment, it may be available

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in the following settings:

- 3 or 5-wire - AAMI standard. Protected against defibrillation.
- 10-wire (optional) - allows up to 12 simultaneous leads. This connector substitutes the standard connector and is not compatible with the 3 or 5-wire cables

7. **Capnography exhaust connector*:** connector used for the removal of the gases collected by capnography.
8. **Capnography connector*:** connector for the capnography sampling line.
9. **NIBP connector*:** connector for direct use with the armband.
10. **SpO₂ connector*:** BCI standard oximetry connector. Adult and child oximetry sensors.
11. **CPR Maestro connector*:** CPR help device connector.
12. **P1 connector:** Invasive pressure connector (channel 1).
13. **P2 connector:** Invasive pressure connector (channel 2).
14. **Defibrillation electrodes (pads) connector:**
 - Multifunctional: adhesive pads for defibrillation, pacemaker and monitoring.
 - Adult/child external pads: accompany the equipment, may be used in adults and children. Cannot be used for pacemaker mode.

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
- Removable battery 1:** the battery can be easily replaced by simply pressing both side tabs one against the other. The battery will unlock and automatically detach itself from the equipment. NOTE: do not remove the battery when the equipment is operating in battery mode. Connect it to an electric current first.
- Removable battery 2:** the battery can be easily replaced by simply pressing both side tabs one against the other. The battery will unlock and automatically detach itself from the equipment. NOTE: do not remove the battery when the equipment is operating in battery mode. Connect it to an electric current first.

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3. **Identification tags:** the identification tags have important information about the product, such as the model, serial number and manufacturer information. This information may be requested if technical assistance is needed. Therefore, do not remove or damage the identification tags.
4. **External DC socket:** for battery connection or external DC source connection in a range of 11 to 16 VDC.
5. **3-pin power connector:** input of 100 to 240 VAC, with central pin for grounding. 5 A fuse (20 mm 20 AG F5A GLASS FUSE).
6. **Grounding and potential equalizer:** potential equalization and general grounding connector.
7. **USB device connector*:** connector for computer connection.
8. **LAN connector*:** connector for wired network connection.
9. **USB host connector*:** connector for pendrive connection.

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 English

Preventive maintenance

Instramed recommends that the DualMax be examined by a qualified technician every 12 months. We recommend that you contact the manufacturer for more information about qualified and trained personnel in your area to perform preventive maintenance.

It is recommended that periodic inspections be performed on the equipment's power supply charger, cables and connectors in order to determine possible isolation or internal conductor ruptures.

Functional tests must be performed at the beginning of every work shift.

Corrective maintenance

If the equipment needs repair, this can only be done by Instramed or its authorized representative, otherwise this Warranty certificate may no longer be valid.

No internal parts are to be fixed by the user.

Removable battery

Even when disconnected (stand-by), the DualMax executes internal routines checking the status of the equipment. Despite of this procedure entailing a low power consumption, the battery charge may be consumed. Therefore, whenever the device has not been connected to an electric current for more than 20 days, it is advisable to execute a full battery charge. If this procedure is not performed, there is a risk of draining the battery and consequently being unable to use the DualMax in its portable configuration (not connected to the electric current). To recharge the battery, connect the monitor to charge the battery, connect the monitor to an AC power source (110 or 220 V outlet) or a DC power source.

There are no restrictions or limitations for using the DualMax while its battery is being recharged by an AC source or DC External source.

Every battery has a determined lifetime, which is the possible quantity of full charge and discharge cycles, without loss of performance (see battery specifications in chapter 8). If the equipment presents a loss in battery performance, please request a new set to Instramed's technical assistance. To request pieces and services, please contact Instramed at +55 (51) 3073 8200.

Battery replacement is recommended every 02 (two) years or when the runtime is less than 01 (one) hour.

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
Removable battery replacement

The battery will automatically detach itself from the equipment as shown in the picture below.

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- 1 - Press the side tabs to unlock the removable battery.
- 2 - Manually remove the battery from the equipment.
- 3 - Correctly position the new battery.
- 4 - Push the new battery until it firmly locks into the cabinet

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Replacing the thermal paper

1 - Pull the lever to open the printer door.

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3 - Place the new bobbin. The bobbin must be positioned as shown in the picture above.

4 - Pull to unroll about 10 cm of paper.

5 - Align paper with the printer door.

6 - Close the printer door. The printer is ready for use

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Disposing of the device

According to the Brazilian environmental legislation, equipments and parts that are no longer in conditions of use should be referred to the manufacturer for the final destination, thus preserving the natural resources and contributing to the conservation of the environment.

For disposal of products from Instramed, contact us by the telephone numbers available on the website www.instramed.com.br or by the e-mail qualidade@instramed.com.br.

Avoid contamination of the environment, humans, or other equipment by making sure to properly decontaminate the equipment before disposing of it.

For countries that follow European Guidelines, refer to 2002/96/CE. For other, countries, refer to local regulations for the proper disposal of trash in your area.

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 English

Instramed Indústria Médico Hospitalar Ltda. warrants the equipment described in this Certificate for 12 (twelve) months, starting from the delivery date. This warranty covers manufacturing or material defects that prevents proper functioning according to the specifications stated herein, as long as the conditions presented in this certificate are respected.

During the warranty period, Instramed Indústria Médico Hospitalar Ltda. or its representative will repair or replace defective parts, at no expense to the equipment's owner.

This warranty will no longer be valid if any damage occurs due to accident, natural disaster, improper connection to a power source, use distinct from that described in the User manual, or irregular working conditions.

Any attempt to violate, adjust or repair this equipment by individuals not authorized by Instramed Indústria Médico Hospitalar Ltda. will automatically invalidate this warranty. This also applies in case of alterations made to this contract, the fiscal receipt, or to the equipment's serial number.

Instramed Indústria Médico Hospitalar Ltda. is not responsible for the improper use of this equipment, by people who are not familiar with its function or the techniques recommend for its proper use.

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ATTENTION: Instramed assumes no responsibility for any damage caused to individuals or property brought by failure to use this product in accordance with the information, recommendations and warnings presented in the user manual, alterations made in the device, attempts of repair not provided by authorized technical assistance centers, operation by unqualified personnel, use of defective device or use of accessories and parts not supplied by the manufacturer.

For information about warranty or technical assistance, please contact Instramed's technical support.

Copyright © 2020 Instramed. The DualMax, Instramed and its respective logos are trademarks of Instramed Indústria Médico Hospitalar Ltda. The internal software of this product is Instramed's intellectual property, being protected under international copyright laws. It is provided exclusively to be used with this present device, identified by the serial number, and may not be, in whole or in part, evaluated, recompiled or altered in any way.

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